





VitaPCR™

Easy. Rapid. Affordable.





Point-of-care Molecular Influenza A/B Assay



2 minutes hands-on time

FEATURES & BENEFITS

✓ Easy

- Intuitive: 2min only hands-on time
- Reliable: no refrigeration needed
- Automated result interpretation

✓ Rapid

- Sample to result in less than 20min
- Streamline your patient management
- Diagnosis and treatment without delay


✓ Affordable

- Same price level of a rapid test
- No trade-off between price & accuracy
- The best diagnosis, for all

SPECS


Sample Collection	Simple Nasal Swab
Sensitivity	2nM
Dimensions	205x165x155mm
Weight	0.9kg

EASY TO USE AND FAST




1. COLLECT

Simple nasal swab sample collection




2. LYSE

Stir the swab in the collection buffer



3. TRANSFER

Use the disposable pipette




4. START

Results delivered in less than 20 minutes


TOTAL PROCEDURE IN LESS THAN 20 MINUTES

IMPACT ON PATIENT MANAGEMENT

Laboratory Batch PCR



VitaPCR





PROVIDING ADVANCED MEDTECH FOR ALL.

Credo Biomedical Pte. Ltd.

1 Fifth Avenue, Guthrie house #02-09 , Singapore 268802

www.credobiomed.com

service@credobiomed.com